

Endotracheal Vacuum Regulator and Vacuum Bottles

Endotracheal Vacuum Regulator are made of Brass Body and suitable to use with Central Medical Vacuum Systems. Available for direct connection into Vacuum Terminal Units with British (BS 5682:1992) DIN (German), French, Italian and American Style Probes.

The Vacuum Regulators can be adjusted at the designed vacuum level.
Vacuum Bottles can be provided in three different models:

- a-) 0.5, 0.85, 1.25 liter direct fitted to Vacuum regulator.
- b-) Wall type 2 liter jar fixed to the wall or rail complete with overflow trap and suction tubes.
- c-) Portable suction trolley having single or double 5 liter suction jar complete, complete with overflow trap, suction tubes and probe.